




The Hoyos Retractor System

*The versatile
abdominal retractor
system that
attaches to your
operating table.*

Proprietary design
facilitates flexibility,
control and
precision

Remora Retractor
designed for
Lipoabdominoplasty


Interlocking
Retractor designed
for Mini-
Lipoabdominoplasty







Dr. Alfredo Hoyos is a world-renowned plastic surgeon from Bogotá, Colombia. He invented Hi-Def and 4D Liposculpture techniques as well as other advanced body contouring techniques. His diverse medical education combined with his background in sculpting have made him a one-of-a-kind plastic surgeon.

Marina Medical is proud to offer superior instrumentation designed and developed with Dr. Hoyos, including the Hoyos Retractor System.



Before 3 months post-op

Remora Retractor for Lipoabdominoplasty



Just like the Remora Fish, who utilizes suction for traction to attach to larger marine animals, the Remora Retractor is designed to atraumatically avoid slippage.

Interlocking Retractors for Mini-Lipoabdominoplasty




The proximal retractor enters through the Cesarean section scars while the distal retractor enters through the umbilicus resulting in a smaller scar and preservation of the umbilicus.

View the Hoyos surgical video and article published in the Aesthetic Surgery Journal at www.marinamedical.com/hoyos

Hoyos Distal Retractor Short
110-603



Hoyos Distal Retractor Long
110-602



Hoyos Proximal Retractor
110-601

